



Name: _____ Department: _____ Student ID: _____

Team No. : _____ Date: _____ Lecturer's Signature: _____

Introduction

1. Standing Waves on A String

Goals

- Understand the transverse wave.
- Understand the superposition principle.
- Observe standing waves on a string.
- Estimate the wave speed on the string.




Figure 1: Standing wave patterns on a string: the first harmonic (top), the second harmonic (middle), and the third harmonic (bottom).




Figure 2: A sonometer system. At certain driving frequencies we observe standing waves.

Theoretical Backgrounds

Consider a transverse vibration along the y axis of a string stretched along the x (longitudinal) axis. For simplicity, we restrict ourselves to the vibration along the y axis and neglect that along the z axis which is actually allowed in nature. We also ignore the vibration along the longitudinal direction.

- (a) A wave is called transverse if the vibration is perpendicular to the longitudinal direction along which the wave propagates.
- (b) A massive string of the linear mass density ρ with a tension τ is a carrier of a propagating wave.
- (c) The speed of the wave propagation along a vibrating string

$$v = \sqrt{\frac{\tau}{\rho}} \tag{1}$$

was discovered in the late 1500s by Vincenzo Galilei, the father of Galileo Galilei.

- (d) A transverse wave is described by a time-dependent field $y(t, x)$, a function of both x and t . Here, $y(t, x)$ is the transverse displacement of a tiny segment of the string at x .
- (e) A wave propagating along the x axis with the propagating speed v is a solution to the wave equation:

$$\frac{1}{v^2} \frac{\partial^2}{\partial t^2} y(t, x) = \frac{\partial^2}{\partial x^2} y(t, x). \tag{2}$$

- (f) For any differentiable function $f(x)$, $y(t, x) = f(x \pm vt)$ is a wave satisfying the wave equation (2). Here, $f(x - vt)$ propagates along $+\hat{i}$ and $f(x + vt)$ along $-\hat{i}$.
- (g) A sinusoidal function

$$y_{\pm}(t, x) = y_m \sin(kx \mp \omega t + \phi) \tag{3}$$

is a wave that propagates along $\pm\hat{i}$ with the speed

$$v = \frac{\omega}{k}.$$

Here, ϕ is a constant phase, ω is the angular frequency and k is the wave number:

$$\omega = \frac{2\pi}{T},$$

$$k = \frac{2\pi}{\lambda},$$

where T and λ are the period and the wavelength, respectively.



- (h) In summary, the wave speed can be expressed as

$$v = \frac{\omega}{k} = \frac{\lambda}{T} = \lambda f,$$

where $f = 1/T$ is the frequency.

- (i) A postulate of the wave is that the resultant wave of any two waves is just the algebraic sum:

$$y(t, x) = y_1(t, x) + y_2(t, x).$$

- (j) We consider the resultant wave of $y_+(t, x)$ and $y_-(t, x)$ that are defined in (3). These waves are of common amplitude and phase except that their directions of propagation are opposite:

$$\begin{aligned} y(t, x) &= y_+(t, x) + y_-(t, x) \\ &= y_m \sin(kx - \omega t + \phi) \\ &+ y_m \sin(kx + \omega t + \phi) \\ &= [2y_m \sin(kx + \phi)] \cos \omega t. \end{aligned} \quad (4)$$

- (k) The resultant wave in (4) does not propagate because the x dependence and the t dependence are factorized (decoupled) completely. Thus the wave does not propagate but stands. Thus we call it a **standing wave**.

- (l) The x dependence $\sin(kx + \phi)$ has zeros periodically at

$$x_{\text{node}} = \frac{n\pi - \phi}{k} = \left(\frac{n}{2} - \frac{\phi}{2\pi}\right) \lambda.$$

We call the points at which the standing wave always vanishes the **nodes**. The distance between the two nearest nodes is $\lambda/2$.

- (m) The x dependence $\sin(kx + \phi)$ has the maximum value 1 periodically at

$$x_{\text{antinode}} = \frac{(n + \frac{1}{2})\pi - \phi}{k} = \left(\frac{n + \frac{1}{2}}{2} - \frac{\phi}{2\pi}\right) \lambda.$$

We call the points at which the standing wave always has the largest vibration the **antinodes**. The distance between the two nearest antinodes is $\lambda/2$.

- (n) For a stretched string of length L with both ends fixed, the two ends must be nodes simultaneously. Therefore, the conditions for the standing waves are

$$L = \frac{\lambda}{2}n, \quad n = 1, 2, \dots$$

- (o) Thus the driver of a sonometer must have the vibrating frequencies satisfying $\lambda_n = 2L/n$ to have standing waves emerge. Thus the wavelengths and frequencies are quantized as

$$f_n = \frac{v}{\lambda_n} = n \frac{v}{2L} = \frac{n}{2L} \sqrt{\frac{\tau}{\rho}},$$

where the harmonic number n is a positive integer and λ_n is the wavelength of the wave at the n th harmonic. We call the lowest frequency (longest wavelength) of $n = 1$ state the **fundamental harmonic** or the **zeroth harmonic**. We call the higher frequency (shorter wavelength) than that of the fundamental harmonic the **overtone**s. For example, $n = 2$ is the first overtone or the first harmonic.


Experimental Procedure

1. E1

- (a) Connect a **PASCO Interface 850** device to a computer.
- (b) Connect a **Detector** to **ANALOG INPUT A** of the interface.




- (c) Select **Sound Sensor** at **Hardware Setup**.







- (d) Connect a **Driver** to **OUTPUT 1** of the interface.
- (e) Apply the AC to the **Driver** by **Signal Generator**.



- (f) Set the **Scope** to measure $V(t)$.



- (g) Set **Fast Monitor Mode** and click **Monitor** to measure the intensity.
- (h) Increase the frequency applied to the **Driver** and find the lowest frequency to maximize the amplitude.




- Note that the resonant frequency of the string is twice the frequency applied to the **Driver**.

- (i) Repeat the step (h) after replacing the density of the string.
- (j) **The harmonic number n (only once during the A8)**: find the second and third lowest frequency to maximize the amplitude without replacing the string.


2. E2

- (a) Connect a **String Vibrator** to **OUTPUT 1** of the interface.
- (b) Connect a **Force Sensor** to **PASPORT 1** of the interface.



- The **PASCO Interface 850** will automatically recognize the **Force Sensor**.


- (c) Set the **Scope** to measure $F(t)$.
- (d) Apply the AC to the **String Vibrator** by **Signal Generator**.
- (e) Set **Fast Monitor Mode** and click **Monitor** to measure the tension of the string.
- (f) Increase the frequency applied to the **String Vibrator** and find the second lowest frequency to maximize the amplitude.



- Skip the frequency when the end of the **String Vibrator** operates as an antinode.
- (g) Initialize a graph of the horizontal and vertical axes with t and $F(t)$, respectively.



- (h) Make linear fits with a fixed slope $m = 0$ to the $F(t)$ data to figure out the tension of the string.



- (i) Replace the tension of the string and find the second lowest frequency to maximize the amplitude.
- (j) Repeat the step (g)-(h).
- (k) **The harmonic number n (only once during the A8)** : find the third and fourth lowest frequency to maximize the amplitude

without replacing the tension.

3. E3

- (a) Connect a **String Vibrator** to **OUTPUT 1** of the interface.
- (b) Apply the AC to the **String Vibrator** by **Signal Generator**.
- (c) Increase the frequency applied to the **String Vibrator** and find the second lowest frequency to maximize the amplitude.
- Skip the frequency when the end of the **String Vibrator** operates as an antinode.
- (d) Replace the length of the string and find the second lowest frequency to maximize the amplitude.
- (e) Repeat the step (d).
- (f) **The harmonic number n (only once during the A8)** : find the third and fourth lowest frequency to maximize the amplitude without replacing the length of the string.



References

[1] **X21.1:** *What is wave?* 

[2] **X21.7:** *Definition of the Wave* 

[3] **X22.1:** *Interference of Waves* 

[4] **X22.2:** *Standing Waves* 

[5] **X22.5:** *Resonance* 